

Specifications of new stepper motor controller RE-C200, RE-C201 for Ring Pump RP-HX & RP-TX & RP-Q-Type

Product name	Low voltage bipolar stepper motor controller
Model number	RE-C200 : Master controller unit RE-C201 : Slave controller unit
Applicable pumps	Ring Pump RP-HX series RP-Q-B series Ring Pump RP-TX series RP-Q-C series RP-HX (D20 motor type - custom spec.)
Drive voltage	RP-HX VM-4.3V(Default VM) / VM-2.6V RP-TX/Q-Type VM-2.6V only
Input voltage	RE-C200 : DC5V 0.5A (USB supply) RE-C201 : DC5V 3.0A (AC-DC Adaptor : Accessory)
Control method	Voltage drive method for bipolar coil.
Excitation method	1/4 Microstep drive method.
Motor speed	1 to 2,000 PPS (4Hz to 8,000Hz) Digital frequency division.
Function	Various controls using PC software 1. SINGLE-MODE 2. PROGRAM-MODE

Connection terminal RE-C200 Size : 100 x 100 x 35mm






Motor connection terminals (CH1) : B/B /A/A

USB : USB mini-Type B terminal (connected to PC with cable bundled.)

Power : LED indicating USB connection

I2C : OUT terminal for connection with slave unit RE-C201

Pump wiring to motor connection terminal



Connection terminal RE-C201 Size : 100 x 160 x 45mm



Address switching SW (A : CH2 to CH5) - (B : CH6 to CH9)

- Motor connection terminal (4CH) : B/B /A/A

I2C : IN terminal for connection with main unit RE-C200 (I2C connection cable)

DC-IN terminal for connection with AC-DC adaptor bundled.

Power : LED indicating power on.

I2C : OUT terminal for connection with slave unit RE-C201

Accessories	RE-C200 :	Type_A-mini_B USB cable USB memory (PC software)
	RE-C201:	AC-DC adaptor I2C connection cable

Motor speed(PPS) and motor rotation by each model

RP-HX series

Motor speed PPS	1	5	10	50	100	200	300	400	500	750	1000		
Motor rotation rpm	0.06	0.30	0.60	3.00	6.00	12.00	18.00	24.00	30.00	45.00	60.00		

RP-TX series

Motor speed PPS	1	5	10	50	100	200	300	400	500	750	1000	1500	2000
Motor rotation rpm	0.02	0.11	0.22	1.10	2.21	4.42	6.63	8.84	11.05	16.57	22.10	33.15	44.20

RP-Q-C Motor Type (1/50)

Motor speed PPS	1	5	10	50	100	200	300	400	500	750	1000		
Motor rotation rpm	0.06	0.30	0.60	3.00	6.00	12.00	18.00	24.00	30.00	45.00	60.00		

RP-Q-B Motor Type (1/150)

Motor speed PPS	1	5	10	50	100	200	300	400	500	750	1000	1100	1200
Motor rotation rpm	0.02	0.10	0.20	1.00	2.00	4.00	6.00	8.00	10.00	15.00	20.00	22.00	24.00

Special notes

- Please refer to the "Motor speed(PPS) & Discharge amount" characteristics for the discharge amount.
(Data document included with the pumps)
- Use the specified control within the range where step-out (motor stop) does not occur.

Single Mode Operation

Single(One-time)simultaneous operation for all connected pumps(up to 9 pumps)

* Connect up to 2 slave units to the main units.

* Each slave unit can control up to 4 pumps.

RE-C200-SGM-V24A-FS1-220119

HELP

<SINGLE-MODE>

▶ **RP-HX** ◀ USB-CONNECT : -OK-

MODEL: RP-HX DIRECTION: CW PPS (Hz):

CH-1 CH-2 CH-3 CH-4 CH-5 **ALL-ON**

CH-6 CH-7 CH-8 CH-9

VOLTAGE: VM-4.3V RPM: 0 ROTATION TO STOP: 1

< ▶▶ FF > <OPERATION> <OPERATION>

START **START** **START**

Until Stop (sec)

Select pump

Click and select from the dropdown options

RP-HX series (Default indication)

RP-TX series

RP-HX-D20 Motor Type (custom model)

RP-Q-B series

RP-Q-C series

Model is selected

Select rotation direction

Click and select from the dropdown options

CW (Clockwise)[Default]

CCW(Counter Clockwise)

Input motor speed

Click and select from the dropdown options

PPS value

1-2,000 pps (Integral value)

Pump and CH selection

CH1 : Main unit (Can be connected with 2 slave units)

Each slave unit can control up to 4 pumps.

Pump and CH selection

Specify ALL-ON/OFF

Voltage selection

RP-HX

Click and select from the dropdown options

VM-4.3V (Default VM) / VM-2.6V

RP-TX/Q-Type

VM-2.6V only

Motor RPM

Indicated automatically calculated motor RPM by designated PPS.

Fast-forward operation

Send liquid or air in the pump tubing by higher speed for faster priming.

Designated PPS motion

START/STOP button to operate at designated PPS settings.

Select number of rotation(S) to automatically stop st.

Click and select from the dropdown options

Select 1 to 15 rotations eg. Select "3" = Pump stops after 3 rotations.

START button to operate at designated number of rotations. Automatically stops.

START motion for ROTATION TO STOP () START Stops automatically.

Operation time indicator (seconds). (Explained in and above)


Approximate time (seconds) from start to stop for Ⓞ.

Program Mode Operation

Programmable operation for all connected pumps (up to 9 pumps)

* Connect up to 2 slave units to the main units

* Each slave unit can control up to 4 pumps.



<Default setting>

Select CHANNEL

Click and select from the dropdown options
Select from CH1 to CH9.

CHANNEL is displayed

*Select the applicable pump connected to the displayed CHANNEL.

Select pump

Click and select from the dropdown options
RP-HX series (Default indication) RP-Q-B series
RP-TX series RP-Q-C series
RP-HX-D20 Motor Type (custom model)

Applicable pump is displayed

Voltage selection

RP-HX
Click and select from the dropdown options
VM-4.3V (Default VM) / VM-2.6V
RP-TX/Q-Type
VM-2.6V only

Fast-Forward (FF) operation

START/STOP for optimal high-speed operation by the applicable pump.
* Select CW(Clockwise)/Default/CCW(Counter Clockwise)
* START to supply by FF with each pumps which connected to each CHANNEL.
Click STOP when solution is delivered to the tip of the tube.

<Program setting>

Select channel Click and select from the dropdown options
Select from CH1 to CH9.

CHANNEL is displayed

* Set the operating conditions for the pump connected to the displayed CHANNEL.

Input motor speed by PPS value

1-2,000 pps (Integrate value)

Select rotation direction

Select CW(Clockwise)/CCW(Counter Clockwise)

ON/OFF time

Set ON (start) time and OFF (stop) time.
It can be set from 00h 00m 01s to 99h 59m 59s.

Setting table (POS)

PPS & RPM Rotation direction ON (start) time & OFF (stop) time are reflected in POS.
Table setting from POS1 to POS3 is possible.

POS-DELETE

Can be cleared after selecting POS1 to POS2. (for re-entering settings)

REPEAT

Specified the number of repeat for setting table POS1 to POS3.
Click and select from the dropdown options
1 to 99 REPEAT time can be selected. (default value is "1")

TIME-TOTAL

Grand total operation time which are setting from POS1 to POS3 of
all CHANNELs ON/OFF of number of repeat.
*Maximum time is 29,700h
Conditions : (99h 59m 59s x3POS) x 99 repeat

Serial Operation

Operate in the order of the selected CHANNEL and STOP when completed.
eg. CH1 setting operation complete CH2 setting operation complete CH3

Parallel operation

Start all the selected channels at the same time and stop when the operation is complete.

Pause

Pause during operation

REFERENCE TIME

Grand total time until stop
*Serial operation Until all channel is completed.
*Parallel operation Until the longest time channel is completed.

Operation mode display

1. Designated channel is shown in black and undesignated channel remains in grey.
2. When designated channel is in operation, channel number is blinking in green.
3. When designated channel is not in operation, channel number is blinking in red.

<Special notes>

START/PAUSE buttons

Buttons will change between "START" and "STOP" and "PAUSE" and "RESTART" after each click.

REFERENCE TIME

Even when the grand total operating time is down to "0" (Zero), it will take 5-15 sec to change from STOP to START, to allow time to complete process completion of all settings.

Program Mode Operation

Programmable operation for all connected pumps (up to 9 pumps)

* Connect up to 2 slave units to the main units

* Each slave unit can control up to 4 pumps.

FILE

RE-C200-PGM-V24A-FS1-220119

FILE HELP

- ALL-RESET
- OPEN-FILE
- SAVE-FILE
- EXIT

<PROGRAM-MODE> USB-CONNECT: -OK-

◀ CHANNEL--1 ▶

CHANNEL: CHANNEL--1 PPS (Hz): 0 <ON>TIME (h:m:s): 0 0 0

REPEAT: 1 DIRECTION: CW CCW <OFF>TIME (h:m:s): 0 0 0

▶	POS.	PPS(Hz)	RPM	DIRECTION	<ON> TIME	<OFF> TIME
▶	1	--	--	--	--:--:--	--:--:--
	2	--	--	--	--:--:--	--:--:--
	3	--	--	--	--:--:--	--:--:--

POS-DELETE TIME TOTAL/ch: -- / 00:00:00 <d/h:m:s>

<SERIAL> START <PARALLEL> START PAUSE

◀ RP-HX ▶

MODEL: RP-HX < >> FF>

VOLTAGE: VM-4.3V FF DIRECTION: CW

CH-1 CH-4 CH-7
CH-2 CH-5 CH-8
CH-3 CH-6 CH-9

<REFERENCE TIME>
Until STOP: -- / 00:00:00 <d/h:m:s>

ALL-RESET
OPEN-FILE
SAVE-FILE
EXIT

All programmed contents are cleared.
Open the programmed and stored file.
Save the programmed file.
Exit.

HELP

RE-C200-PGM-V24A-FS1-220119

FILE HELP

- MANUAL
- VERSION

<PROGRAM-MODE> USB-CONNECT: -OK-

◀ CHANNEL--1 ▶

CHANNEL: CHANNEL--1 PPS (Hz): 0 <ON>TIME (h:m:s): 0 0 0

REPEAT: 1 DIRECTION: CW CCW <OFF>TIME (h:m:s): 0 0 0

▶	POS.	PPS(Hz)	RPM	DIRECTION	<ON> TIME	<OFF> TIME
▶	1	--	--	--	--:--:--	--:--:--
	2	--	--	--	--:--:--	--:~:~:~
	3	--	--	--	--:~:~:~	--:~:~:~

POS-DELETE TIME TOTAL/ch: -- / 00:00:00 <d/h:m:s>

<SERIAL> START <PARALLEL> START PAUSE

◀ RP-HX ▶

MODEL: RP-HX < >> FF>

VOLTAGE: VM-4.3V FF DIRECTION: CW

CH-1 CH-4 CH-7
CH-2 CH-5 CH-8
CH-3 CH-6 CH-9

<REFERENCE TIME>
Until STOP: -- / 00:00:00 <d/h:m:s>

MANUAL
VERSION

Manual of software in Japanese and in English
Version of the program