


## RE-C602 Controller Specifications for Ring Pump RP-Q-F Bi-polar Stepper Motors

<b>Product Name</b>	Low voltage bi-polar stepper motor controller
<b>Model Number</b>	RE-C602 : controller unit
<b>Applicable Pumps</b>	Ring Pump RP-Q-F
<b>Drive Voltage</b>	RP-Q-FL Series: VM=3.3 V (RP-Q / F1.5S-2P6A-DC3.3VS) RP-Q-FH Series: VM=4.3 V (RP-Q / F1.5S-3P6A-DC4.3VS)
<b>Input Voltage</b>	RE-C602 : DC5V, 1.0 A (AC-DC adaptor of AC100-240V 50/60Hz is included as accessory)
<b>Control Method</b>	Voltage drive method for bi-polar coil
<b>Excitation Method</b>	1/4 microstep drive method
<b>Motor Speed</b>	1 - 2,000 pps (4 - 8,000 Hz) Digital frequency division
<b>Functions</b>	Various controls using PC software <b>SINGL-MODE</b> <b>Program-Mode</b>
<b>Connection Terminals:</b>	RE-C602 (Size: 100 x 100 x 34 mm)



- USB Connection Terminal (USB A-microB: Included as accessory)
- Pump Connection Terminals-1
- Pump Connection Terminals-2
- DC-IN Connection Terminal (AC Adaptor: Included as accessory)
- Start Switch
- Start-LED (BLUE)
- Pause-Switch
- Pause-LED (GREEN)
- Power-LED (RED)



- Accessories Included:** RE-C602 :
- Type A-micro B USB Cable
  - USB Memory Stick (PC Software)
  - AC-DC Adapter


## Motor Speed (PPS) and Motor Rotation for Each Series

RP-Q-F Series (1/21.06 gear ratio motor specification)												VM = 4.3V	
Motor Speed PPS	1	5	10	50	100	200	300	500	800	1000	1100	1200	1400
Motor Rotation rpm	0.14	0.71	1.42	7.12	14.25	28.49	42.74	71.23	113.96	142.45	156.75	171.00	199.43

### Special Notes:

- a) Discharge Rate: Please input "Discharge Rate per Rotation ( $\mu\text{L}/\text{rev}$ )".  
(Refer to data included with the pumps)
- b) Control the motors within the range to prevent motor stepping-out (or motor stopping or vibrating).

# SINGLE




## <Single Settings>

- ① Pump Selection
  - Click the "v" (drop down) and select the pump.
  - RP-Q-FL
  - RP-Q-FH
  - Operating Voltage: VM value is selected automatically.
  - RP-Q-FL VM = 3.3V (fixed)
  - RP-Q-FH VM = 4.3V (fixed)
- ② DIRECTION
  - rotating direction. (Clockwise CW or Counter-Clockwise CCW)
- Pump Speed
  - 1 ~ 2 0 0 0 PPS (Limit:1400)
- FREE
  - Click FREE to supply media (solution) for pump
  - Click STOP.

## <PROGRAM Settings>

- Input Resolution "uL/rev"
  - Input the value from "Discharge Volume per Rotation (uL/rev)" that comes attached with the pump.
- ENTER
  - Click ENTER to fix the "Resolution" settings.
- ROTATION
  - 1 ~ 3 0
- Flow Rate(uL)
  - ${}^1\text{uL/rev}_1 * {}^1\text{ROTATION}_1 = \text{Flow Rate (uL)}$
- START
  - START button while connected to the PC (computer) with USB.
  - \* Auto Stop
- EXIT
  - Disconnect from PC (computer).

# PROGRAM



## <Initial Settings>

- |                            |   |
|----------------------------|---|
| Pump Selection             | Click the "▼" (drop down) and select the pump.<br>RP-Q-FL<br>RP-Q-FH  |
| VOLTAGE                    | Operating Voltage: VM value is selected automatically.<br>RP-Q-FL VM = 3.3V (fixed)<br>RP-Q-FH VM = 4.3V (fixed)  |
| Fast Forward Operation - 1 | Click the "▼" (drop down) and rotating direction. (Clockwise CW or Counter-Clockwise CCW)   |
| Fast Forward Operation - 2 | START/STOP the optimal high speed pump operation to supply the media (solution).<br>* Click START to supply media (solution) for each pump on each channel. Once the media (solution) flows to the tip of the discharge tube, click STOP. |

## <PROGRAM Settings>

- |                            |   |
|----------------------------|---|
| Input Resolution "µL/rev"  | Input the value from "Discharge Volume per Rotation (µL/rev)" that comes attached with the pump.  |
| Specify Flow Rate "µL/min" | Input desired flow rate (µL/min).   |
| ENTER                      | Click ENTER to fix the "Resolution" and "Flow Rate" settings. The optimal "PPS" and "RPM" values will be automatically displayed.   |
| Select DIRECTION           | Select between Clockwise (CW) or Counter-Clockwise (CCW)  |
| Operation Time             | Set operation (ON) time<br>* Setting range between 00h 00m 01s to 99h 59m 59s   |
| Stop Time                  | Set stop (OFF) time<br>* Setting range between 00h 00m 01s to 99h 59m 59s<br>(3)  |
| POS-CLEAR                  | Select between POS 1 to POS 5 ( ) and click POS-CLEAR to clear setting for that POS. (To reset settings)  |
| RESET                      | RESET all settings in POS 1 to POS 5.   |
| Specify POS                | Click the "▼" (drop down) and select POS.<br>Allows up to 5 different program condition settings in POS 1 to POS 5.<br>Program is reflected on the POS by repeating to .<br>( Flow Rate, ENTER, DIRECTION, Operation Time, Stop Time) |
| REPEAT                     | Input number of repeat operation for specified program table (POS 1 to POS 5).<br>Initial setting is "1". Can input value between 1 to 99.  |
| WAIT TIME (h:m)            | If needed, specify the wait time before the START time for each channel. (Input h:m)  |

## <WRITE Settings>

- |       |   |
|-------|---|
| WRITE | Click WRITE to input the PROGRAM settings to controller's ROM (memory). (During WRITE, it will be in WAIT mode) |
|-------|---|

**<Operation Settings>**

**"Controller Switch - START"**

Use the controller's switch to START / STOP / PAUSE without the PC.


**"PC-START"**

- START
- PAUSE
- EXIT

START / STOP button while connected to the PC (computer) with USB.  
 \* Each click will change between START STOP.  
 PAUSE button while connected to the PC (computer) with USB.  
 \* Each click will change between PAUSE RESTART PAUSE.  
 Disconnect from PC (computer).

**NB. When controlling from PC side, "START" and "Pause" of main unit SW must be set to <OFF> position to prevent possible unstable operation.**

**SAVE-FILE / OPEN FILE**



SAVE-FILE

Click on SAVE-FILE to save the PROGRAM setting (Type in file name)

OPEN-FILE

Click on OPEN-FILE to open saved PROGRAM setting.

- \* Once the PROGRAM setting is opened, click WRITE to save it in the computer ROM.
- \* Click all the CH that needs to operate (letters turn from grey to black), and click CH-SET to save it in the computer ROM (memory).