Ring Pump®

RP-G II Series

Features

- Compact box shape makes it possible to pile up plural G II pumps.
- We have 2G II series which can flow two different liquids by a single pump.


Examples of application

Automatic washing machine Electrolytic water conditioner Applications for detergent supply etc.

Main specifications


Flow rate (mL/min) *	10 12 17 20 25 40 45			
Tubing	Silicone/ PharMed BPT/ Tygon A-60-F			
Motor/Power source	DC geared motor DC12V、24V			
Weight (for reference)	Approx. 90g			

[Outline view]



[Mounting plate view]


[Flow direction]



Model Lineup

Model Lineup	Voltage (DC/V)	Flow rate (mL/min)	Tubing	
			Туре	ID x OD(mm)
RP-G II 02 S -10Z-DC12V	12	12	Silicone	2 x 4
RP-G II 02F-10Z-DC24V	24	10	PharMed BPT	2 x 4
RP-G II 02F-10Z-DC12V	12	10	PharMed BPT	2 x 4
RP-G II 02 S -25Z-DC24V	24	25	Silicone	2 x 4
RP-G II 02F-25Z-DC24V	24	25	PharMed BPT	2 x 4
RP-G II 2.5S-20Z-DC12V	12	20	Silicone	2.5 x 4.5
RP-G II 2.5S-40Z-DC24V	24	40	Silicone	2.5 x 4.5
RP-G II 2.8N-20Z-DC12V	12	20	Tygon A-60-F	2.8 x 4.8
RP-G II 2.8N-40Z-DC24V	24	40	Tygon A-60-F	2.8 x 4.8

^{*} PharMed BPT and Tygon A-60-F are products of Saint-Gobain SA and formerly described as PharMed and Norprene in our brochure.



Variable by material or ID of tubing.